




Press Release FG08-10

Fused Silica Capillary Tubing from Fiberguide Industries

Polyimide coating can withstand temperatures up to 360°C

Stirling, N.J. (August 17, 2010) – Delivering reliable operation in demanding conditions, Fiberguide Industries introduces its new line of fused silica capillary tubing. The new line of tubing has proven effective in a number of scientific applications including protein characterization, small molecule and genetic analysis, proteomics, DNA sequencing and fragment analysis.



Click on this low resolution image to download a full resolution file

Fiberguide's capillary tubing is manufactured from high purity synthetic quartz with a protective coating of polyimide. Applied to the outer surface of the tubing, polyimide has an upper temperature limit of 360°C. A wide variety of standard I.D.'s are available, from 2 μ to 700 μ and O.D.'s from 90 μ to 850 μ . Non-standard sizes can also be manufactured by request.

For additional information on these and other Fiberguide products and services, including new data sheets to assist with fiber selection, contact Fiberguide toll free at 877-490-7803, email info@fiberguide.com or visit www.fiberguide.com.

Fiberguide Industries, Inc. manufactures a comprehensive line of standard and custom high optical transmission fibers, OEM assemblies and ultra precision arrays. The company is a part of the global Halma technology group, which also includes photonics specialists Ocean Optics (www.oceanoptics.com), a leader in spectroscopy and optical thin films; and Labsphere (www.labsphere.com), a recognized expert in the light metrology field. FDA registered as a Contract Manufacturer and Custom Device Manufacturer, Fiberguide's corporate and optical fiber manufacturing facilities are located in Stirling, New Jersey, with a manufacturing/assembly facility in Caldwell, Idaho.

###

To download a full resolution file of the image, please click on the image above or visit: http://halmapr.com/fiberguide/FG_CapillaryTubing_EMAIL.jpg

This Press Release and an archive of all Fiberguide news can also be found on the Fiberguide News Blog at: <http://halmapr.com/news/fiberguide>

Enquiries to:

Konrad Lisi
Fiberguide Industries
1 Bay Street
Stirling, NJ 07980 USA
Tel: +1 908.647.6601
Fax: +1 908.647.8464
E-mail: info@fiberguide.com
Website: www.fiberguide.com

Press Contact:

Rob Dietrich
Halma Holdings PR Services
11500 Northlake Dr., Suite 306
Cincinnati, OH 45249
Tel: +1 513.898.8007
Fax: +1 513.898.8008
E-mail: rob.dietrich@halma.com
Website: www.halmapr.com